



Laparoscopic Anti Fog Solution

Clarify your vision



gri

alleset[®]

DELIVERING PRODUCT SOLUTIONS[™]

www.gri-alleset.com



Used by surgeons to reduce fog on laparoscopes, our Anti Fog Solution offers clear results. With a separate 6-gram bottle and sponge in a convenient kit, the sponge can be easily re-saturated as needed throughout the surgical case.

The sponge is a non-abrasive foam that is safe for use on today's delicate scopes and includes a x-ray detectible element for added safety.

Feature	Benefits
<i>6-Gram Bottle with Dropper</i>	<i>Plenty of solution to easily moisten the sponge for even the longest surgical procedures</i>
<i>Adhesive Backed Sponge</i>	<i>Allows for convenient and secure placement anywhere on the sterile field</i>
<i>Non-Abrasive Sponge</i>	<i>Reduces possibility of damaging delicate scope lenses</i>
<i>X-Ray Detectible Element in Foam Pad</i>	<i>Assists in x-ray detection if lost in the surgical site</i>
<i>Does not contain natural rubber latex or phthalates</i>	<i>Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds</i>

Description	REF	Case Qty.
Anti-Fog Solution 6ml Bottle Non-sterile*	37-001	2000 / case
Sponge Adhesive Pad for Anti-Fog Solution with X-Ray Element Non-sterile	37-002	2000 / case
Anti-Fog Solution Kit 6ml Bottle with Sponge Sterile (gamma)	37-005-S	200 / case

* Non-sterile bottle has been pre-sterilized via gamma irradiation. Liquid inside of bottle is validated as sterile and once assembled and re-packaged, it may be EO sterilized.

Corporate Office

GRI-Alleset
4142 Industry Way
Flowery Branch, Georgia 30542 USA
Main: +1.678.866.0550
Fax: +1.678.866.0551
Email: info@gri-alleset.com
Web: www.gri-alleset.com



DELIVERING PRODUCT SOLUTIONS™
www.gri-alleset.com