



allesorb™
SOLIDIFIER

Solidifying fluids. Solidifying safety.



gri™

alleset®

DELIVERING PRODUCT SOLUTIONS™

www.gri-alleset.com



Waste fluid management is an ongoing concern in any healthcare facility. Anything that can help prevent hazardous spills resulting from leaking, cracked, or unsealed suction canisters will be a welcome product.

GRI-Alleaset Allesorb Solidifier is a high-performance superabsorbent polymer that, when added to a suction canister, quickly gels the contained liquids to help prevent potentially hazardous spills. Allesorb Solidifier sets the standard in fast and effective waste fluid control.

Feature	Benefits
<i>Yorker cap</i>	<i>“Funnel” style cap can be used when only suction ports are available – otherwise, just unscrew the end cap and use the pour spout for faster dispensing</i>
<i>Two Common Sized Bottles</i>	<i>Sized according to the suction canisters you use – no need to measure, just pour in the entire contents for a pre-measured volume.</i>
<i>Convenient Dissolvable Pouch</i>	<i>For smaller volumes of waste fluid, just drop-in the pouch. No need to tear and pour: the pouch will dissolve in the liquid.</i>
<i>No Shaking Required</i>	<i>Allesorb will gel liquid without the need to shake or stir the waste container, helping to reduce the chances of spills or splashes.</i>
<i>Non-Sandy Solid Gel</i>	<i>Allesorb gels in larger clumps, rather than a sandy texture. If a suction canister breaks, the contents of a more solid mass are easier to clean-up</i>
<i>Does not contain natural rubber latex or phthalates</i>	<i>Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds.</i>



The dissolvable pouch is a simple and effective alternative to bottles. Simply drop the pouch through the suction canister pour spout, and watch as the paper dissolves, simultaneously releasing the solidifier powder and efficiently gelling the liquid.

Description	REF Non-Sterile	Case Qty.
<i>Solidifies 1000cc, Dissolvable Pouch (33gm)</i>	<i>70-S1000</i>	<i>100 / Case</i>
<i>Solidifies 1500cc, Bottle (1.6 oz / 45 gm)</i>	<i>70-S1500</i>	<i>96 / Case</i>
<i>Solidifies 3000cc, Bottle (3.2 oz / 91 gm)</i>	<i>70-S3000</i>	<i>96 / Case</i>

Corporate Office

GRI-Alleaset
 4142 Industry Way
 Flowery Branch, Georgia 30542 USA
 Main: +1.678.866.0550
 Fax: +1.678.866.0551
 Email: info@gri-alleaset.com
 Web: www.gri-alleaset.com



DELIVERING PRODUCT SOLUTIONS™
www.gri-alleaset.com